



Heine Alpha+ 3.5v Fiber Optic Otoscope

- Distal Fiber Optic illumination with XHL Xenon Halogen Technology.
- Homogenous, reflex-free illumination of ear canal and tympanum.
- Optical glass viewing window with 3x magnification.
- Scratch-resistant, high-resolution image.
- Swivelling viewing window, permanently attached.
- Ideal for instrumentation, cannot get mislaid.
- Insufflation port for pneumatic testing.
- High-impact polycarbonate housing with chrome-plated cone for durability.

D-02.80.114

\$445.00


WelchAllyn® PocketScope Otoscope

- 2.5V Otoscope also has built in throat illuminator - two in one instrument! Features halogen illumination and fiber optics for a cool light with no reflections or obstructions, plus a rectangular viewing lens for instrumentation under magnification.
- 2.5V PocketScope handle has a smooth rheostat to control light intensity and features a handy pocket clip. Key Components
- PocketScope™ Otoscope (Head Only)
- 2.5v PocketScope™ Handle (without battery)
- 2 - "AA" Battery * 2.5v Halogen Lamp

\$163.99

22820



Heine mini 3000 Otoscope

Modern pocket otoscope in a unique, compact design. Can only be used with the mini 3000 handle system.

- Enhanced XHL Xenon Halogen Technology with 100% more light. Very bright, concentrated white light for perfect illumination.
- Viewing window with 3x magnification and optimized casing surface for razor-sharp images and minimal reflection.
- Swivelling viewing window built into instrument. Useful for instrumentation, cannot be mislaid.
- Attachment clip with integrated on/off switch. Secure. Switches off automatically when returned to the pocket.
- 2-piece handle and head. Easy to maintain, flexible.
- High-quality handle: Chrome-finish upper section/ refined plastic. Shockproof, sturdy, non-slip.
- Replaceable batteries. Size AA.

D-01.70.220

\$135.00


Heine Mini 3000 Fiber Optic Otoscope Set

- Mini 3000 F.O. Otoscope
- 1 set = 4 reusable tips
- 5 each of 2.5 and 4 mm dia. AllSpec disposable tips
- Hard case

D-851.10.021

\$257.00

WelchAllyn® 3.5v MacroView Diagnostic Otoscope Set



25070M

\$299.00

- Twice the field of view – you can see almost the whole of the tympanic membrane without panning
- 30% greater magnification than a traditional otoscope – from 2.5x to 4.2x
- Better clarity – for improved definition of landmarks
- Adjustable focusing – suits different lengths of ear canal and compensates for far-sighted eyes
- Choice of specula – disposable or reusable available in a range of sizes
- Insufflation port – for pneumatic otoscopy