



Hillrom™

**LET US
CHANGE
YOUR VIEW**

Explore Welch Allyn Physical Exam Sets and Systems featuring our new, LED wide-view PanOptic Plus Ophthalmoscopes and MacroView Plus Otoscopes

WALL SYSTEMS AND TRANSFORMERS

			OPHTHALMOSCOPES			OTOSCOPES & ACCESSORIES				
	Wall Board	777 Wall Transformer	PanOptic Plus with Quick Eye Alignment	PanOptic Basic	Coaxial	MacroView Plus, Digitally-Capable	MacroView Basic	iExaminer SmartBracket for MacroView Plus	Ear Spec Dispenser	Wall Aneroid
			118-3-US	118-2-US	11720-L	238-3	238-2	11840-IEX-US		
777 Integrated Wall Systems	777-PM3WCX-US	✓	✓	✓		✓		✓	✓	
	777-PM3WAS-US	✓	✓	✓		✓		✓	✓	✓
	777-PM3WAX-US	✓	✓	✓		✓		✓	✓	✓
	777-PM3WXS-US	✓	✓	✓		✓		✓	✓	
	777-PM3WXX-US	✓	✓	✓		✓		✓	✓	
	777-PM2WCX-US	✓	✓		✓		✓		✓	
	777-PM2WAS-US	✓	✓		✓		✓		✓	✓
	777-PM2WAX-US	✓	✓		✓		✓		✓	✓
	777-PM2WXS-US	✓	✓		✓		✓		✓	
	777-PM2WXX-US	✓	✓		✓		✓		✓	
Wall Transformers	777-SM2WXX	✓	✓		✓		✓		✓	
	777-SM2WAS	✓	✓		✓		✓		✓	✓
	777-SM2WAX	✓	✓		✓		✓		✓	✓
	777-PM3XXX-US		✓	✓		✓		✓		
	777-PM2XXX-US		✓		✓		✓			
	777-SM2XXX		✓		✓		✓			

*Choosing these models will ensure the necessary mounting bracket ships with your wall board.

Welch Allyn Connex Spot Monitor and/or Spot Vital Signs Device must be ordered separately.

**5-year warranty is only applicable to the PanOptic Plus Ophthalmoscope and MacroView Plus Otoscope heads.

DIAGNOSTIC SETS

	OPHTHALMOSCOPES			OTOSCOPES & ACCESSORIES				POWER HANDLES (sets include one power handle)		
	PanOptic Plus with Quick Eye Alignment	PanOptic Basic	Coaxial	MacroView Plus, Digitally-Capable	MacroView Basic	Traditional	iExaminer SmartBracket for MacroView Plus	Lithium-Ion Plus with built-in USB	Lithium-Ion	Nickel Cadmium
	118-3-US	118-2-US	11720-L	238-3	238-2	250-2	11840-IEX-US			
71-PM3LXE-US	✓			✓			✓	✓		
71-XM3LXE				✓			✓	✓		
71-PX3LXERT-US	✓		✓					✓		
71-PX3LXE-US	✓							✓		
71-PM3LXES-US	✓			✓			✓	✓		
71-PM2LXU-US		✓			✓				✓	
71-SM2LXU			✓		✓				✓	
71-SS2CXX			✓			✓				✓
71-SM2CXX			✓		✓					✓

**5-year warranty is only applicable to the PanOptic Plus Ophthalmoscope and MacroView Plus Otoscope heads.

VITAL SIGNS

SureTemp Plus Thermometer	Connex Spot Monitor/4400 Vital Signs Device*	5-Year Warranty**
	✓	✓
✓		✓
✓		✓
✓		✓
✓	✓	✓
✓		✓
✓		✓
		✓

Welch Allyn Green Series 777 Integrated Wall Systems and Transformers can help you stay efficient during your examinations. Select from wide-view and standard ophthalmoscopes and otoscopes, digital or manual blood pressure devices and integrated or stand-alone oral/axillary/rectal thermometry with a **Welch Allyn SureTemp Plus Thermometer**.



#777-PM3WCX-US

Model	Hard Case	Soft Case	5-Year Warranty**
	✓		✓
	✓		✓
	✓		✓
	✓		✓
		✓	✓
	✓		✓
	✓		✓
	✓		✓

Welch Allyn Diagnostic Sets are the perfect solution to keep your exam tools secure for whenever you need them. Try our wide-view optics and fast charge power supplies built into our Plus models.



#71-PM3LXE-US

DESK SETS

All Welch Allyn Desk Sets shown below ship with a Universal Desk Charger and power supply.

Welch Allyn Desk Sets are compact and perfect for smaller offices to support eye and ear examinations. Try our wide-view optics our Plus models.




#71-PM3LDE-US

	OPHTHALMOSCOPES			OTOSCOPES & ACCESSORIES			POWER HANDLES (sets include two of the same handle)		
	PanOptic Plus with Quick Eye Alignment	PanOptic Basic	Coaxial	MacroView Plus, Digitally-Capable	MacroView Basic	iExaminer SmartBracket for MacroView Plus	Lithium-Ion	Nickel Cadmium	5-Year Warranty**
	118-3-US	118-2-US	11720-L	238-3	238-2	11840-1EX-US			
71-PM3LDE-US	✓			✓		✓	✓		✓
71-PM2LDX-US		✓			✓		✓		
71-SM2LDX			✓		✓		✓		
71-SM2NDX			✓		✓			✓	
71-XXXLDE							✓		

**5-year warranty is only applicable to the PanOptic Plus Ophthalmoscope and MacroView Plus Otoscope heads.

Welch Allyn Ophthalmoscopes



	11720-L	118-2-US	118-3-US
Model	Traditional	PanOptic Basic	PanOptic Plus
Warranty	1 year	1 year	5 year
Larger field of view	—	20X	20X
Designed with DuraShock	—	Yes	Yes
Quick Eye alignment technology	—	—	Yes

Welch Allyn Oscopes



	250-2	238-2	238-3
Model	Traditional	MacroView Basic	MacroView Plus
Warranty	1 year	1 year	5 year
Larger field of view	—	3X	3X
Designed with DuraShock	—	Yes	Yes
Digitally-capable for image capture	—	—	Yes

Connect with your Baxter representative today to discuss the best Welch Allyn Physical Exam & Diagnostic solution for your clinical needs.



hillrom.com

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products. Baxter, Connex, DuraShock, Hillrom, iExaminer, MacroView, PanOptic, Quick Eye, SmartBracket, SureTemp and Welch Allyn are trademarks of Baxter International Inc. or its subsidiaries. ALL RIGHTS RESERVED. APR226002 R3 EN-US 22-APR-2022