# **!-!**

# Hillrom...





# SURETEMP PLUS 690 SETTING THE STANDARD OF CARE

Advance your temperature routine with the **SureTemp** Plus 690 Thermometer, featuring technology that sets the standard for speed and accuracy.



An easy-to-read LCD display with flashing icons that clearly communicate in any language, and an intuitive user interface and ergonomic shape.

# 2 SAFE

Patient safety is improved thanks to an interchangeable oral/axillary and rectal probe wells that help reduce the risk of cross-contamination.

## **3** FAST AND ACCURATE

In addition to improved accuracy, the **SureTemp** Plus 690 is one of the fastest thermometers on the market, providing approximately 4-6 second oral readings, 10-13 second pediatric axillary and rectal readings and 12-15 second adult axillary readings.

### **4** VERSATILE

To top it all off, you can now use the **SureTemp** Plus 690 with all your patients, from infants to adults—whether you need to take oral, axillary or rectal temperatures with an easy-to-use, intuitive operation.



# **OPTIONAL WALL MOUNT**

The **SureTemp** Plus 690 with optional wall mount has storage space for 50 additional probe covers.





#### SETTING THE STANDARD OF CARE EVEN HIGHER

#### **TECHNICAL SPECIFICATIONS**

Calibration Accuracy	±0.2 °F (±0.1 °C) (Monitor Mode)
Temperature Predict Time	Oral: Approximately 4-6 seconds Adult Axillary (18 years and older): Approximately 12 - 15 seconds Pediatric Axillary (17 years and younger): Approximately 10 - 13 seconds Rectal: Approximately 10 - 13 seconds
Unit Dimensions	8.46" x 3.18" x 2.43" (215 mm x 81 mm x 62 mm)
Weight	12 oz (357 g)
Ambient Temperature Range	50.0 °F (10 °C) to 104 °F (40.0 °C)
Patient Temperature Range	80.0 °F (26.7 °C) to 110 °F (43.3 °C)
Power Source	(3 ) 1.5 Vdc AA batteries
Battery Life	Approximately 6,000 temperatures

#### PRODUCT AND ACCESSORY ORDERING INFORMATION

#### **THERMOMETERS**

**01690-200:** SureTemp Plus 690 Electronic Thermometer, Oral Probe and Oral Probe Well with 4 ft cord

**01690-201**: SureTemp Plus 690 Electronic Thermometer, Rectal Probe and Rectal Probe Well with 4 ft cord

**01690-300:** SureTemp Plus 690 Electronic Thermometer, One per room, Wall Mount, Oral Probe and Oral Probe Well with 9 ft cord (must be locked to the wall)

**01690-301:** SureTemp Plus 690 Electronic Thermometer, One Per Room, Wall Mount, Rectal Probe and Rectal Probe Well with 9 ft cord (must be locked to the wall)

**01690-400:** SureTemp Plus 690 Electronic Thermometer, Wall Mount, Oral Probe and Oral Probe Well with 4 ft cord

01690-700: SureTemp Plus 690 Electronic

#### **ACCESSORIES**

21326-0000: Wall Holder

21330-0000: Holder for Spare Probe and Well 06138-000: Calibration Key Assembly, 690/692 02892-000: Probe and Well Kit, 4ft, Rectal 02892-100: Probe and Well Kit, 9ft, Rectal 02893-000: Probe and Well Kit, 4ft, Oral 02893-100: Probe and Well Kit, 9ft, Oral

**02891-0000:** Probe Well, Blue **02891-1000:** Probe Well, Red

#### **REPLACEMENT PROBE COVERS**

**05031-101:** Case of 1.000 Probe Covers

For more information, please contact your local distributor or Hillrom sales representative at 1-800-535-6663.

#### hillrom.com

Baxter, Hillrom, SureTemp, and Welch Allyn are trademarks of Baxter International, Inc. or its subsidiaries.